

Contact Ensures Reliable Connections for Implantables

The Bal Conn contact, which consists of a canted-coil spring retained in a metal housing, ensures uniform, reliable mechanical and electrical contact between the battery and lead in implantable devices. The compact, innovative design of the Bal Conn enhances performance and reduces the connection envelope for the critical devices used in cardiac healthcare management and neurostimulation therapies. Bal Conn designs are available that support the emerging IS4 cardiac standard.

Source URL (retrieved on 02/01/2015 - 11:45am):

<http://www.mdtmag.com/product-releases/2008/11/contact-ensures-reliable-connections-implantables>