

Peristaltic Filling System

The Multi Filling System elevates peristaltic technology into a smooth, seamless process in which the sterile product comes into contact with only the medically graded silicone tubing and filling needle. When the process is completed, the tubing can be easily cleaned, or treated as a disposable, to eliminate contamination and facilitate the clean-up validation process.

Source URL (retrieved on 01/25/2015 - 11:53pm):

http://www.mdtmag.com/product-releases/2004/02/peristaltic-filling-system?qt-recent_content=0&qt-video_of_the_day=0