

TPR Tubing (1513)

Suprene tubing is produced under Good Manufacturing Practices (GMP) using first-quality polypropylene-based TPR. It offers excellent low-temperature flexibility, is lighter in weight than comparable sizes of rubber tubing, and provides good resistance to tearing and fatigue from repeated flexing. It also offers outstanding compression characteristics and durability, making it an ideal fit for use with peristaltic pumps. To meet various requirements, Suprene is available in both FDA and industrial grades.

Source URL (retrieved on 01/25/2015 - 10:20pm):

<http://www.mdtmag.com/product-releases/2004/12/tpr-tubing-1513>